

Varianta a typ svítidla	Objednací číslo	Světelný výkon / 1 m (lx)	Teplota chromatičnosti (K)	Index barev - CRI (Ra)	Průměr světél. pole (cm)	Hloubka osvětlení (cm)	Délka ramene (cm)	Hmotnost (kg)
Stropní								
TRIANGO 100-1 C	ND15.884.000	100 000	4 300	> 96	18	80	152	15,9
TRIANGO 100-1 C Fokus	ND15.886.000	100 000	4 300	> 96	18 / 23 / 28	80	152	15,9
TRIANGO 100-1 C Endo	ND15.972.000	100 000	4 300 / Endo mode	> 96	18	80	152	15,9
TRIANGO 100-3 C	ND15.885.000	100 000	3 700 / 4 300 / 4 700	> 96	18	80	152	15,9
TRIANGO 100-3 C Fokus	ND15.887.000	100 000	3 700 / 4 300 / 4 700	> 96	18 / 23 / 28	80	152	15,9
Nástěnné								
TRIANGO 100-1 W	ND15.888.000	100 000	4 300	> 96	18	80	152	15,3
TRIANGO 100-1 W Fokus	ND15.890.000	100 000	4 300	> 96	18 / 23 / 28	80	152	15,3
TRIANGO 100-1 W Endo	ND15.973.000	100 000	4 300 / Endo mode	> 96	18	80	152	15,3
TRIANGO 100-3 W	ND15.889.000	100 000	3 700 / 4 300 / 4 700	> 96	18	80	152	15,3
TRIANGO 100-3 W Fokus	ND15.891.000	100 000	3 700 / 4 300 / 4 700	> 96	18 / 23 / 28	80	152	15,3
Pojízdné								
TRIANGO 100-1 F	ND15.892.000	100 000	4 300	> 96	18	80	72	20
TRIANGO 100-1 F Fokus	ND15.894.000	100 000	4 300	> 96	18 / 23 / 28	80	72	20
TRIANGO 100-1 F Endo	ND15.974.000	100 000	4 300 / Endo mode	> 96	18	80	72	20
TRIANGO 100-3 F	ND15.893.000	100 000	3 700 / 4 300 / 4 700	> 96	18	80	72	20
TRIANGO 100-3 F Fokus	ND15.895.000	100 000	3 700 / 4 300 / 4 700	> 96	18 / 23 / 28	80	72	20
TRIANGO 100-1 F Battery	ND16.222.000	100 000	4 300	> 96	18	80	72	32
TRIANGO 100-3 F Fokus Battery	ND16.223.000	100 000	3 700 / 4 300 / 4 700	> 96	18 / 23 / 28	80	72	32