

HALUX N30-1 F



TECHNICAL DATA

FITTED WITH	1 x LLED module
LUMINOUS POWER	30 000 Lux / 0,5 m *
ILLUMINANCE FIELD 0.5 m	ø 18 cm*
COLOR RENDERING INDEX	CRI > 93
COLOR TEMPERATURE	4 300 k*
OPERATION EQUIPMENT	Electronic converter
VOLTAGE	100-240 V; 50/60 Hz
CORD LENGTH	2 m
LUMINAIRE BODY	Polycarbonate (PC) & acrylonitrile butadiene styrene copolymer (ABS)
WEIGHT (NET)	4,4 kg
POWER CONSUMPTION	8 W
USAGE	Push button on light head
CLASS OF PROTECTION	II
PROTECTION RATING	IP 20
STANDARDS MET	EN 60601-1, EN 60601-2-41
SPECIAL FEATURES	Continuously dimmable, maintenance-free
MOUNTING	Roller stand
ACCESSORIES	D16.064.000-00808869

EXAMINATION LIGHTS




LIGHT DISTRIBUTION



TECHNICAL DRAWING

