



fitted with	LED module - Energy efficiency category A+
Illuminancefield @ 50 cm	Ø 18 cm
illuminance @ 50 cm	30 000 lx
color rendering index (CRI)	> 93
colour temperatur	4400 K
work equipment	electronic plug-in transformer
mains lead	1.7 m AUS, Euro, UK, US
luminaire body	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
weight (net)	1.6 kg
Power consumption	8 W
usage	keypad at luminaire head
balance of articulation joints	friction
class of protection	II
Accessory	roll support, Rail fastening with removal guard and brake function, rail fastening with removal protection, table clamp with removal protection, universal fastening with removal protection, extension arm with removal protection, wall bracket with removal protection
norms	EN 60601-1, EN 60601-2-41
special features	continuously dimmable, maintenance free
Fastening	Fastening pin (Ø 16mm)
order number	D16156000 - pure white

light distribution

