



<b>fitted with</b>	1 x LED module - Energy efficiency category A+
<b>Illuminancefield @ 50 cm</b>	Ø 18 cm
<b>illuminance @ 50 cm</b>	50 000 lx
<b>color rendering index (CRI)</b>	> 93
<b>colour temperatur</b>	4400 K, 3800 K, 3300 K
<b>work equipment</b>	electronic converter
<b>mains lead</b>	1.8 m Euro
<b>luminaire body</b>	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
<b>weight (net)</b>	1.3 kg
<b>Power consumption</b>	18 W
<b>usage</b>	keypad at luminaire head
<b>balance of articulation joints</b>	friction
<b>class of protection</b>	II
<b>Accessory</b>	roll support, Rail fastening with removal guard and brake function, rail fastening with removal protection, table clamp with removal protection, universal fastening with removal protection, extension arm with removal protection, wall bracket with removal protection
<b>norms</b>	EN 60601-1, EN 60601-2-41
<b>special features</b>	3 color shades, 5 dimming levels, maintenance free
<b>Fastening</b>	Fastening pin (Ø 16mm)
<b>order number</b>	D16157000 - pure white

## light distribution

