

Minor surgical light

TRIANGO 100-1 C DUO D16201000-00811242

Waldmann **W**

ENGINEER OF LIGHT.



Derungs 
LIGHT FOR PEOPLE

fitted with	LED - Energy efficiency category A+
illuminancefield @ 100 cm	Ø 18 cm (each)
illuminance @ 100 cm	100 000 lx (each)
color rendering index (CRI)	> 95 (each)
colour temperatur	4300 K (each)
work equipment	electronic converter
mains lead	direct power supply (2x)
luminaire body	polyamide (PA)
weight (net)	30.0 kg
Power consumption	80 W
usage	keypad at luminaire head
balance of articulation joints	spring
class of protection	I
Accessory	sterilizable handle, adapter for suspended ceiling
norms	EN 60601-1, EN 60601-2-41
special features	4 dimming levels, Diamond optics, High-quality arm system, sterilisable grip, maintenance free
Fastening	ceiling mount, when mounting ceiling tube could be cut to length 25-120 cm
System of protection	IP 20 (43)
order number	D16201000- pure white

light distribution

