

TRIANGO 100-1 C

TECHNICAL DATA



FITTED WITH	1 x LED module
LUMINOUS POWER	100 000 lux / 1.0 m *
COLOR RENDERING INDEX	Ra > 95
OPERATION EQUIPMENT	Electronic converter
VOLTAGE	100 – 240 V; 50/60 Hz
CORD LENGTH	Direct power supply
LUMINAIRE BODY	Polyamide (PA)
WEIGHT (NET)	15.9 kg
POWER CONSUMPTION	40 W
USAGE	Touch panel on cardanic system
TECHNOLOGY	Single-pole secondary
BALANCE OF ARTICULATION JOINTS	Spring (weight compensation)
CLASS OF PROTECTION	I
PROTECTION RATING	IP 43 (luminaire head)
STANDARDS MET	EN 60601-1, EN 60601-2-41
SPECIAL FEATURES	Diamond optics, sterilizable handle, maintenance-free, high quality arm system, dimmable
MOUNTING	Ceiling mount, ceiling tube can be cut on site (length: 25-120 cm)
ACCESSORIES	Sterilizable handle, intermediate ceiling adapter

TREATMENT LIGHT




	illumination area / 1.0 m		color temperature		dimmable	Order No.
	fixed value 18 cm	changeable value 18 / 23 / 28 cm	fixed value 4300 K	changeable value 3700 / 4300 / 4700 K		
TRIANGO 100-1 C	x				x	D15.884.000-00802695
TRIANGO 100-3 C	x			x	x	D15.885.000-00803113
TRIANGO Fokus 100-1		x	x		x	D15.886.000-00803685
TRIANGO Fokus 100-3 C		x		x	x	D15.887.000-00803686
TRIANGO Endo 100-1 C	x		x		x Endo mode	D15.972.000-00803704

TECHNICAL DRAWING



valid for Triango 100 with fixed illumination area



measurements in centimeters * -10% / +20% tolerance