



TRIANGO 60 C



|                                       |  |
|---------------------------------------|--|
| <b>Fitted with</b>                    | LED - Energy efficiency category A+<br>ø 18 cm                                   |
| <b>Illuminance 100 cm</b>             | 60 000 lx  |
| <b>Colour rendering index (Ra)</b>    | > 93   |
| <b>Work equipment</b>                 | Electronic converter   |
| <b>Voltage</b>                        | 100-240 V; 50/60 Hz  |
| <b>Mains lead</b>                     | Direct power supply  |
| <b>Luminaire body</b>                 | Polyamid (PA)  |
| <b>Weight (net)</b>                   | 13 kg  |
| <b>Power consumption</b>              | < 30 W   |
| <b>Usage</b>                          | Keypad at cardanic   |
| <b>Balance of articulation joints</b> | Spring   |
| <b>Class of protection</b>            | I  |
| <b>Accessory</b>                      | Sterilizable handle, handle cover,<br>adapter for suspended ceiling              |
| <b>Norms</b>                          | EN 60601-1, EN 60601-2-41  |
| <b>Special features</b>               | Stepless dimmable,<br>diamond optics, sterilizable handle                        |
| <b>Fastening</b>                      | Ceiling mount, when mounting<br>ceiling tube could be cut to<br>length 25-120 cm |
| <b>System of protection</b>           | IP 43 (luminaire head)   |
| <b>Medical device class</b>           | I according to MDR   |

|                         | lighting field /<br>1,0 m<br><br>fixed value<br>18 cm | lighting field /<br>1,0 m<br><br>changeable<br>value<br>18 - 25 cm | Colour<br>temperature<br><br>fixed value<br>4700K | Colour<br>temperature<br><br>changeable<br>value<br>3700K / 4300K<br>/ 4700K | Dimmable | Order No.            |
|-------------------------|---|--|---|--|----------|----------------------|
| TRIANGO 60-1 C          | X   |  | X   |  | X        | D16.250.000-00813741 |
| TRIANGO 60-1 C US       | X   |  | X   |  | X        | D16.253.000          |
| TRIANGO 60-3 C          | X   |  |   | X  | X        | D16.251.000          |
| TRIANGO 60-3 C US       | X   |  |   | X  | X        | D16.254.000          |
| TRIANGO Fokus 60-1 C    |   | X  | X   |  | X        | D16.252.000          |
| TRIANGO Fokus 60-1 C US |   | X  | X   |  | X        | D16.255.000          |

Light distribution



Technical drawing

