



## TRIANGO 60-1 F



<b>Fitted with</b>	LED - Energy efficiency category A+
<b>Illuminance field 100 cm</b>	ø 18 cm
<b>Illuminance 100 cm</b>	60 000 lx
<b>Colour rendering index (Ra)</b>	> 93
<b>Work equipment</b>	Electronic converter
<b>Voltage</b>	100-240 V; 50/60 Hz
<b>Mains lead</b>	Direct power supply
<b>Luminaire body</b>	Polyamid (PA)
<b>Weight (net)</b>	19 kg
<b>Power consumption</b>	< 30 W
<b>Usage</b>	Keypad at cardanic
<b>Balance of articulation joints</b>	Spring
<b>Class of protection</b>	I
<b>Accessory</b>	Sterilizable handle, handle cover
<b>Norms</b>	EN 60601-1, EN 60601-2-41
<b>Special features</b>	Stepless dimmable, diamond optics, sterilizable handle
<b>Fastening</b>	Roller stand
<b>System of protection</b>	IP 43 (luminaire head)
<b>Medical device class</b>	I according to MDR

	Lighting field / 1,0 m fixed value 18 cm	lighting field / 1,0 m changeable value 18 - 25 cm	Colour temperature fixed value 4700K	Colour temperature changeable value 3700K / 4300K / 4700K	Dimmable	Order No.
TRIANGO 60-1 F CH	X		X		X	D16.262.000-00813744
TRIANGO 60-1 F Schuko	X		X		X	D16.259.000-00813743
TRIANGO 60-1 F UK	X		X		X	D16.271.000-00813747
TRIANGO 60-1 F AUS	X		X		X	D16.268.000
TRIANGO 60-1 F US	X		X		X	D16.265.000

### Light distribution



### Technical drawing

