





Fitted with	LED - Energy efficiency category A+
Illuminance field 100 cm	∅ 18 cm
Illuminance 100 cm	60 000 lx
Colour rendering index (Ra)	> 93
Work equipment	Electronic converter
Voltage	100-240 V; 50/60 Hz
Mains lead	Direct power supply
Luminaire body	Polyamid (PA)
Weight (net)	12,5 kg
Power consumption	< 30 W
Usage	Keypad at cardanic
Balance of articulation joints	Spring
Class of protection	I
Accessory	Sterilizable handle, handle cover
Norms	EN 60601-1, EN 60601-2-41
Special features	Stepless dimmable, diamond optics, sterilizable handle
Fastening	Wall mount
System of protection	IP 43 (luminaire head)
Medical device class	I according to MDR

	lighting field / 1,0 m	lighting field / 1,0 m	Colour temperature	Colour temperature	Dimmable	Order No.
	fixed value 18 cm	changeable value 18 - 25 cm	fixed value 4700K	changeable value 3700K / 4300K / 4700K		
TRIANGO 60-1 W	X		X		X	D16.256.000-00813742
TRIANGO 60-1 W US	X		X		X	D16.274.000
TRIANGO 60-3 W	X			X	X	D16.257.000
TRIANGO 60-3 W US	X			X	X	D16.275.000
TRIANGO Fokus 60-1 W		X	X		X	D16.258.000
TRIANGO Fokus 60-1 C US		X	X		X	D16.276.000

Light distribution



Technical drawing

